

FOR IMMEDIATE RELEASE – August 25, 2008

Contact: Bob Luton, President & CTO
Bio Image Systems, Inc.
721 Bloomfield
Jackson, MI 49203 USA

**Advanced Quantifier and Intelligent Quantifier®
Available for Apple Mac OS X**

Bio Image Systems, Inc. announces *Advance Quantifier* and *Intelligent Quantifier®* are available for Apple Macintosh computers running OS X with Intel (Mactel). You can download a free, fully functional demo version from our website, www.bioimage.net. The download file includes instructions for installing on Mac. The installation requires installation of VMware Fusion, or a similar program such as Parallels Desktop or Apple Boot Camp.

Advanced Quantifier is a complete software system for analyzing scanned images of 1-D wet gels, dry gels, and membranes: lane definition, band detection, Mw/ IEF standards, and dilution series (concentration standards). You can match both within and between samples: to all lanes, to a selected (reference) lane, to selected bands, or to a reference pattern. Matching output includes a similarity matrix and a dendrogram.

Intelligent Quantifier® analyzes images from 1-D wet gels, dry gels, and membranes. It has a complete one-dimensional analysis package including lane definition, band detection, Mw/ IEF standards, and dilution series (concentration standards). The software generates band and lane data analysis reports. *Intelligent Quantifier®* does not have lane matching capabilities.

About Bio Image Systems, Inc.

Bio Image Systems, Inc. has been developing and selling high-quality image analysis software and systems since 1982. Bio Image products are available for Windows PCs and OS X Apple Macintosh.

For more information, contact Bob Luton at Bio Image Systems, Inc.

Phone: +1.734.604.2623

Fax: +1.517.788.8869

email: info@bioimage.net

Internet: <http://www.bioimage.net>