

FOR IMMEDIATE RELEASE – June 6, 2008

Contact: Bob Luton, President & CTO
Bio Image Systems, Inc.
721 Bloomfield
Jackson, MI 49203 USA

Intelligent Quantifier-11 Version 1.2 Adds Interface to Excel® Macros

Bio Image Systems, Inc. announces Version 1.2 of *Intelligent Quantifier-11 for Windows 2000 and XP*, with an interface to Excel® macros. An easy-to-use template combines lane and band data with a Visual Basic processor to load any Excel® spreadsheet. The processor can report up to 20 appropriate loading errors, which are displayed both on the monitor and entered into the record's audit trail.

Intelligent Quantifier-11, a software system for analyzing scanned images of 1-D electrophoresis separations, is fully compliant with 21CFR Part 11. The software includes fully integrated logins and passwords; record signings (acquire image, save changes, and reports); Analyst, Supervisor, and Administrator access levels; secure audit trails of all operations and a system audit trail.

The *Intelligent Quantifier-11* software analyzes images from wet gels, dry gels, and membranes. It has a complete one-dimensional analysis package including lane definition, band detection, Mw/IEF standards, and dilution series (concentration standards). The software generates band and lane data analysis reports and comes with an IQ/OQ. Customized reports and complete hardware and software systems are also available.

About Bio Image Systems, Inc.

Bio Image Systems, Inc. has been developing and selling high-quality image analysis software and systems since 1982. Bio Image products are available for Windows PCs and Sun workstations.

For more information, contact Bob Luton at Bio Image Systems, Inc.

Phone: +1.734.604.2623

Fax: +1.517.788.8869

email: info@bioimage.net

Internet: <http://www.bioimage.net>