



Bio Image Systems, Inc.

721 Bloomfield
Jackson, MI 49203-3959 USA

Phone: +1.734.604.2623
FAX: +1.517.788.8869
Email: info@bioimage.net
<http://www.bioimage.net>

ACS 5000 Specifications

Light Source:

Wavelength	400-750 nm
Color	User-selected Gray, Red, Green, or Blue
Mode	Transmissive and Reflective Absorbance
Range	0.05-3.77 (transmissive), 0-2.0 RD (reflective)
Bit Depth	16 bit (65536 gray levels)
Power	100-240V, 50-60 Hz

Scan Area:

Spatial Resolution	User-selected 100, 150, 300, or 600 dpi
Sample Area (cm)	20 x 25 (transmissive) 30 x 20 (reflective)
Dimensions (cm)	30 W x 16 H x 50 D