

## FOR IMMEDIATE RELEASE – May 8, 2006

Contact: Bob Luton, President & CTO  
Bio Image Systems, Inc.  
721 Bloomfield  
Jackson, MI 49203 USA

### **Intelligent Quantifier-11 Version 1.1 1-D Electrophoresis Analysis Software Provides 21CFR Part 11 Compliance**

Bio Image Systems, Inc. announces Version 1.1 of *Intelligent Quantifier-11*, the first commercially available image analysis system fully compliant with 21CFR Part 11. The *Intelligent Quantifier-11* software analyzes images from wet gels, dry gels, and membranes.



New features of Version 1.1 include

- The QImaging Retiga® firewire camera joins the Roper Scientific Sensys® camera and Microtek i900® flatbed scanner as a directly supported acquisition device. Device settings are stored with the metadata of each record.
- A streamlined calibration process.
- Multiple lanes per report page, with customizable headers and footers for both pages and lanes.
- 5 new items in the header or footer. *Record Signing Comments* is a quick way to add a report-specific comment to each page.

The software includes fully integrated logins and passwords; record signings (acquire image, save changes, and reports); Analyst, Supervisor, and Administrator access levels; secure audit trails of all operations and a system audit trail.

It has a complete one-dimensional analysis package including lane definition, band detection, Mw/ IEF standards, and dilution series (concentration standards). The

software generates band and lane data analysis reports with fully customizable reports in saveable templates.

There is an optional 75-page IQ/OQ. Customized reports and complete hardware and software systems are also available. The software runs on Windows 2000 and Windows XP.

### **About Bio Image Systems, Inc.**

Bio Image Systems, Inc. is a sales and support company for **Bio Image** software. For over 20 years the name Bio Image has meant powerful, yet easy-to-use image analysis software to thousands in the pharmaceutical and biotechnology industries.

Today Bio Image continues to be a leader with *Intelligent Quantifier-11*, the first commercially available image analysis program for electrophoresis samples which is **fully compliant with 21CFR Part 11**.

At Bio Image Systems, Inc. we understand our customers' needs - whether in QA/QC, methods development, or research - and can provide not only the best imaging solution possible, but also the high quality support desired. We are proud of our products, and proud of our reputation in the industry.

For more information, contact Bob Luton at Bio Image Systems, Inc.

**Phone:** +1.734.604.2623

**Fax:** +1.517.788.8869

**email:** [info@bioimage.net](mailto:info@bioimage.net)

**Internet:** <http://www.bioimage.net>